#### NH2-750F **TOYOPEARL®** SALT TOLERANT

#447


# INTRODUCTION

**EXCHANGE RESIN** 

Ion Exchange Chromatography (IEC) is one of the most frequently used chromatographic modes for the separation and purification of biomolecules. Compared with other chromatographic modes, modern ion exchange media offer high dynamic binding capacities and a straightforward method development. IEC is used at all stages and scales of purification of therapeutic proteins: from laboratory scale purification to industrial scale downstream processing. Today, modern ion exchange resins offer extremely high binding capacities. The interest is now shifting towards salttolerant ion exchange media that enable capturing out of a biological feedstock at physiological conditions or direct processing of target fractions without dilution.

Tosoh Bioscience has developed a new salt tolerant anion exchange resin, TOYOPEARL NH2-750F, offering a high binding capacity across a range of pH values and conductivities. A TOYOPEARL HW-75 polymeric bead has been functionalized with a primary amine-containing ligand resulting in a resin with increased salt tolerance and selectivity different than that of quaternary amine anion exchange resins currently available. This new resin is ideal for process scale applications from the capture of proteins from biological feedstock (mammalian cell culture, plasma, bacterial feedstock, etc.) without dilution to the intermediate or final purification of monoclonal antibodies (mAbs) where aggre-

COMPARISON OF SALT TOLERANCE OF AEX RESINS



Resins: various; Column size: 5 mm ID × 5 cm L; Mobile phase: A: 20 mmol/L Tris-HCl, pH 8.0, B: mobile phase A + 2.0 mol/L NaCl; Gradient: 0 -100% B (120 min); Flow rate: 300 cm/h (1.0 mL/min); Detection: UV @ 280 nm; Temperature: ambient; Sample: BSA (1.0 g/L)

gates and other impurities are removed from the target of interest. TOYOPEARL NH2-750F is capable of post-protein A removal of aggregates in both, flow-through and bind-elute modes.

### **HIGHLIGHTS**

- ➤ High protein binding capacity at elevated ion strength
- ➤ Unique selectivity, differing from conventional AEX media
- Binds targets at pH close to isoelectric point
- ➤ Aggregate removal in flow-through or bind/elute mode

TOYOPEARL NH2-750F offers static binding capacities approaching 70 g/L for bovine serum albumin across a range of pH values and conductivities. Increased salt tolerance of TOYOPEARL NH2-750F as compared to other TOYOPEARL anion exchange resins can be seen in Figure 1. While BSA is starting to elute at 0.14 mol/L NaCl for most conventional ion exchangers, the BSA peak begins to elute from the TOYOPEARL NH2-750F column at a concentration of approximately 1 mol/L NaCl.

Retention can be affected by mobile phase pH (Figure 2) without greatly changing the selectivity of the resin. BSA binding ocuurs even at pH values similar to the isoelectric point, indicating a multimodal binding mechanism of the resin. This allows for a large design space in which to develop a separation method.

# INFLUENCE OF pH ON BSA ELUTION




Figure 2

Resins: TOYOPEARL NH2-750F; Column size: 5 mm ID x 5 cm L; Mobile phase: A: 20 mmol/L N-methyl piperazine, pH 4.5 and 5.0; 20 mmol/L Bis-Tris, pH 6.0; 20 mmol/L Tris-HCl, pH 7.5 B: mobile phase A + 2.0 mol/L NaCl; Gradient: 0 -100% B (120 min); Flow rate: 300 cm/hr (1.0 mL/min); Detection: UV @ 280 nm; Temperature: ambient; Sample: BSA (1 mg, pl 4.7-4.9)






TOYOPEARL NH<sub>2</sub>-750F is effective at removing aggregates from mAbs, as demonstrated in Figure 3. SEC analysis of the peaks (data not shown) shows that fraction 1 contains pure monomer. High molecular weight aggregates are completely removed from the mAb peak.


TOYOPEARL NH<sub>2</sub>-750F offers good caustic stability and exhibits excellent pressure-flow characteristics (Figure 4).

## SEPARATION OF AGGREGATES FROM IgG<sub>1</sub> MONOMER



Resin: TOYOPEARL NH<sub>2</sub>-750F; Column size: 5 mm ID  $\times$  5 cm L; Mobile phase: A: 20 mmol/L Tris-HCl, pH 8.0; B: mobile phase A + 1.0 mol/L NaCl; Gradient: 0 -100% B (60 min); Flow rate: 300 cm/h (1.0 mL/min); Detection: UV @ 280 nm; Temperature: ambient; Sample: mAb (lgG<sub>1</sub>) 0.5 g/L

#### PRESSURE - FLOW CHARACTERISTICS OF TOYOPEARL NH2-750F



Resin: TOYOPEARL NH2-750F; Column size: 4.4 cm ID  $\times$  29 cm L; Mobile phase: 0.1 mol/L NaCl; Flow rate: multiple



### **Ordering Information**

#### Part-No Description Resin volume Pore size Particle size **TOYOPEARL** 0023438 TOYOPEARL NH2-750F 100 mL > 100 nm 45 µm TOYOPEARL NH2-750F 250 mL 0023439 > 100 nm 45 µm 0023440 TOYOPEARL NH2-750F 1 L > 100 nm 45 µm 0023441 TOYOPEARL NH2-750F 5 L > 100 nm 45 µm 0023442 TOYOPEARL NH2-750F 50 L > 100 nm 45 µm **ToyoScreen** 0023443 ToyoScreen NH2-750F 1 mL x 6 > 100 nm 45 µm 0023444 ToyoScreen NH2-750F 5 mL x 6 > 100 nm 45 µm MiniChrom 0045108 MiniChrom TOYOPEARL NH2-750F

## Headquarters

Tramstraat 15 5611 CM Eindhoven

**RoboColumns** 

0045021

0045022

**TOYOPEARL NH2-750F** 

5611 CM Eindhoven T +31 (0) 40 257 39 72 F +31 (0) 40 251 47 58

# Sales and Service

Apolloweg 2B 8239 DA Lelystad T +31 (0) 320 87 00 18

Box 3 1831 Diegem T +32 (0) 2 721 92 11



5 mL

200 µl x 8

600 µl x 8

**TOSOH BIOSCIENCE** 

> 100 nm

> 100 nm

> 100 nm

45 µm

45 µm

45 µm

With courtesy of

Max-Planck-Strasse 4 D-47475 Kamp-Lintfort T +49 (0) 28 42 9280 799 F +49 (0) 28 42 9732 638

ToyoScreen RoboColumn NH2-750F

ToyoScreen RoboColumn NH2-750F

Cedar Court, Grove Park Business Est. White Waltham, Maidenhead Berks, SL6 3LW T +44 (0) 16 288 220 48 F +44 (0) 70 394 006 78

info@go-jsb.com www.go-jsb.com

