

Water Management in Purge & Trap/GC

One of the most crucial features of purge and trap/GC instrumentation is the water management system. The *CDS Analytical PeakMaster* is a purge and trap system specifically designed to minimize water management problems.

There are two components to the system. The first is a WETrap that consists of 3/8 inch stainless steel tubing between the purge bottle and the absorbent trap inside a temperature-controlled insulator. During the purge time, the WETrap remains at ambient temperature so that water condenses inside the tubing, and doesn't enter the absorbent trap. The advantages of the WETrap are shown in Figure 1. 5 ml of water was purged with helium at 40 ml/min for a purge time of 11 minutes, with the sample at 85 C. The amount of water that entered the GC column with no WETrap present was measured and compared to the amount that entered the GC column with a WETrap. 98.9% of the water that was carried by the helium sweep gas was removed with the WETrap.

The second component of the water management system consists of the trap dry function, during which the absorbent trap is vented at ambient temperature to allow the water to be flushed through. The trap is then rapidly heated during the desorb function to sweep the analytes onto the GC column. As shown in Figure 2, after as little as a 1 minute trap dry period, only 5 % of the water that entered a Tenax absorbent trap remained to enter the GC column. This means that with the combined use of the WETrap and trap dry func-

Figure 1





Figure 2



tions, 99.95% of the water is removed before the sample enters the GC column.

With the use of the PeakMaster WETrap and trap dry functions, water problems are virtually eliminated during applications requiring a Tenax trap. However, the silica gel required by the EPA to trap gases during EPA 502.2 and EPA 601


analyses poses a particular problem. Figure 3 shows the effect of water retention on the EPA 502.2 gases separation. The first chromatogram was obtained at the beginning of the day after

long bake period for the trap, and the second chromatogram shows the same sample under the same conditions after a 24-hour period of analyses using an autosampler. Note that peak 4, bromomethane, is no longer detected in the second chromatogram. The deterioration of the separation is due to a gradual build up of water in the system.

To ensure that water accumulation doesn't cause a problem when using an absorbent trap containing silica gel, we recommend that the trap dry function last 5 minutes during each analysis, and that the

trap bake function last 5 to 10 minutes. Once a day, the trap should be baked at 200 C for 30 minutes to an hour. This will ensure that no gradual accumulation of water occurs.

Figure 3



© courtesy of CDS Analytical Inc.

Headquarters

JSB International
Tramstraat 15
5611 CM Eindhoven
T +31 (0) 40 251 47 53
F +31 (0) 40 251 47 58

Zoex Europe
Tramstraat 15
5611 CM Eindhoven
T +31 (0) 40 257 39 72
F +31 (0) 40 251 47 58

Sales and Service

Netherlands
Apolloweg 2B
8239 DA Lelystad
T +31 (0) 320 87 00 18
F +31 (0) 320 87 00 19

Belgium
Grensstraat 7
Box 3 1831 Diegem
T +32 (0) 2 721 92 11
F +32 (0) 2 720 76 22

Germany
Max-Planck-Strasse 4
D-47475 Kamp-Lintfort
T +49 (0) 28 42 9280 799
F +49 (0) 28 42 9732 638

UK & Ireland
Cedar Court,
Grove Park Business Est.
White Waltham, Maidenhead
Berks, SL6 3LW
T +44 (0) 16 288 220 48
F +44 (0) 70 394 006 78

info@go-jsb.com
www.go-jsb.com

