


A Practical Applications Guide for Analytical Pyrolysis – GC/MS


Forensic



CDS Analytical www.cdsanalytical.com

Powder-Coat Paints

Placing a suspect at the scene of a crime is crucial when solving crimes. This is achievable through the identity of trace evidence. Evidence such as paints, fibers, toner, ink and cosmetics can easily be characterized. In this example, powder-coat paints are differentiated from each other. Epoxy (top pyrogram) is identified by bisphenol A, and isopropenyl phenol, while benzoic acid is indicative of polyester (bottom pyrogram). A blend or hybrid (center pyrogram) has pyrolysis products of each.




Typical powder coatings. Epoxy (top), Polyester (bottom) and Hybrid (center).

Pyroprobe SettingSet-point:750°C for 15 sValve Oven:325°CTransfer Line:325°CGC/MS Setting		
Valve Oven:325°CTransfer Line:325°CGC/MS SettingColumn:35% PhenylCarrier:Helium, 50:1 splitInjector:325°COven:40°C for 2 min 10°C/min to 300°C	Pyroprobe Setting	
Transfer Line:325°CGC/MS SettingColumn:35% PhenylCarrier:Helium, 50:1 splitInjector:325°COven:40°C for 2 min 10°C/min to 300°C	Set-point:	750°C for 15 s
GC/MS Setting Column: 35% Phenyl Carrier: Helium, 50:1 split Injector: 325°C Oven: 40°C for 2 min 10°C/min to 300°C	Valve Oven:	325°C
Column:35% PhenylCarrier:Helium, 50:1 splitInjector:325°COven:40°C for 2 min 10°C/min to 300°C	Transfer Line:	325°C
Column:35% PhenylCarrier:Helium, 50:1 splitInjector:325°COven:40°C for 2 min 10°C/min to 300°C		
Carrier:Helium, 50:1 splitInjector:325°COven:40°C for 2 min 10°C/min to 300°C	GC/MS Setting	
Injector:325°COven:40°C for 2 min10°C/min to 300°C	Column:	35% Phenyl
Oven: 40°C for 2 min 10°C/min to 300°C	Carrier:	Helium, 50:1 split
10°C/min to 300°C	Injector:	325°C
	Oven:	40°C for 2 min
Scan Range: 35-550 amu		10°C/min to 300°C
	Scan Range:	35-550 amu



Toner Inks

Toner materials used in laser printers and photocopiers are frequently a combination of organic polymers or copolymers and inorganic materials like metal oxides used as pigments. Toners can be analyzed printed onto paper, important for questioned document investigations. The following pyrograms are toners printed on paper. Pyrolysis products of the toner, peaks 1,2,4, and 5, are acrylate monomer and styrene oligomers are easily recognizable from the paper itself.




- 1. Styrene monomer
- 2. Butyl acrylate monomer
- 3. Levoglucosan (from paper)
- 4. Styrene dimer
- 5. Styrene trimer



Textile Fibers

Fibers transferred to a crime scene can be strong evidence. Often, one fiber is plenty of sample to identify its polymer type. Fiber blends are like other polymer blends or mixtures when analyzed; the individual polymers pyrolyze independently, and the resulting pyrogram contains pyrolysates from each of the constituents. In the following example, trace fibers both contained in a fabric listed as 30% cotton. When pyrolyzed, Fiber 1 has only pyrolysis products for PET polyester, but Fiber 2 has peaks from cellulose in cotton (specifically furfural and levoglucosan) as well as the benzoic and terephthalic acid products from PET.



Cotton/polyester fiber blends. Fiber 1 is pure polyethylene terephthalate (PET). Fiber 2 is a blend of PET and cotton.

Pyroprobe Setting	
Set-point:	750°C for 15 s
Valve Oven:	300°C
Transfer Line:	325°C
GC/MS Setting	
Column:	5% phenyl
Carrier:	Helium, 75:1 split
Injector:	300°C
Oven:	40°C for 2 min
	10°C/min to 325°C
Scan Range:	35-550 amu

Peak # Compound

- 1 Vinyl benzoate
- 2 Benzoic acid
- 3 Divinyl terephthalate
- 4 Terephthalic acid monovinyl ester
- 5 Furfural
- 6 Levoglucosan
- 7 Acrylonitrile
- 8 Hexene dinitrile



CDS Analytical Headquarters:

465 Limestone Road P.O. Box 277 Oxford, PA 19363-0277 Tel: 800.541.6593 610.932.3636 © CDS Analytical LLC. 2018 All Rights Reserved AP-12 11/08 5K Printed in USA