

Agilent Sulfur Chemiluminescence Detector and Nitrogen Chemiluminescence Detector

Specification Guide

Agilent 8355 SCD

 MDL^1 Typical <0.5 pg(S)/sec

Typical selectivity g S/g C^2 >2 × 10⁷ Linearity³ >10⁴

Repeatability <2% RSD 2 hours <5% RSD 24 hours

Ozone supply gas Dry oxygen (Ultra Zero grade)

Oxidizer Ultra Zero grade air Hydrogen Ultra Zero grade

Analog output 0-1 V, 0-10 V (Standalone version only)

- 1 MDL: Burner temperature 800 °C, 8 mL/min lower hydrogen, 38 mL/m upper hydrogen, 50 mL/min air, tert-butyl disulfide in SCD checkout sample (5190-7003) as the test compound, 30 m × 0.32 mm, 1 μ m DB1 (123-1033), 50 °C for 3 minutes, 25 °C/min to 160 °C, hold 2 minutes, 1 μ L splitless injection, fully stabilized burner.
- ² Selectivity: Defined as the sensitivity of S over the sensitivity of a selected hydrocarbon. Operating parameters same as MDL. Isooctane (sample solvent) and tert-butyl disulfide in the SCD checkout sample are used as the testing compounds.
- 3 Linearity: Operating parameters same as MDL except 12 mL/min lower hydrogen and 42 mL/m upper hydrogen; 80 °C for 1 minute, 25 °C/min to 160 °C, hold 2 minutes, COC inlet; Test compound: tert-butyl disulfide in isooctane.



Agilent 8255 NCD

 MDL^1 Typical <3 pg(N)/sec in both N and

nitrosamine modes

Typical selectivity g N/g C^2 > 2 × 10⁷ in N mode (selectivity in

nitrosamine mode is matrix-dependent)

Linearity³ >10⁴

Repeatability <1.5% RSD 8 hours

(~ the same in N and nitrosamine) < 2% RSD 18 hours (~3% RSD in nitrosamine mode over 21 hours)

Ozone supply gas Dry oxygen (Ultra Zero grade)

Oxidizer Ultra Zero grade oxygen

Hydrogen Ultra Zero grade

Analog output 0-1 V, 0-10 V (Standalone version only)

- 1 **MDL**: Burner temperature 900 °C, 3 mL/min hydrogen, 8 mL/m oxygen, 3-methylindole in NCD checkout sample (5190-7002) as the test compound, 30 m × 0.32 mm, 0.25 μm HP-5 (19019J-413), 50 °C for 3 minutes, 25 °C/min to 250 °C, hold 2 minutes, 1 μL splitless injection, fully stabilized burner. (For nitrosamine mode, dipropylnitrosamine in methanol is a good test compound, COC inlet, 0.5 μL injection, a WAX column is preferred to match the polarity of nitrosamines).
- 2 Selectivity: Defined as the sensitivity of N over the sensitivity of a selected hydrocarbon. Operating parameters same as MDL. Isooctane (sample solvent) and 3-methylindole in the NCD checkout sample are used as the testing compounds.
- 3 **Linearity:** Operating parameters same as MDL except oven temperatures (80 °C for 1 minute, 25 °C/min to 180 °C, hold 1 minute, nitrobenzene in isooctane as test compound.

Physical Specifications

Power requirements

8255/8355 Detector and Pump 120/220-240 V

50/60 Hz 1.200 VA

Dimensions

Detector Height: 41.0 cm (16.1 in)

Width: 27.0 cm (10.6 in)
Depth: 51.1 cm (20.1 in)

8355 SCD weight 22 kg (49 lbs) 8255 NCD weight 24 kg (52 lbs)

Burner Height: 22.1 cm (8.7 in)

Weight: 0.7 kg (1.5 lbs)

Vacuum pump Height: 26.1 cm (10.3 in)

Width: 15.8 cm (6.2 in) Depth: 43.0 cm (16.9 in) Weight: 25 kg (55 lbs)

Environmental conditions

Installation category II
Pollution degree 2

Ambient temperature 50–104 °F (10–40 °C) Relative humidity 80% at 87.5 °F (31 °C)

50% at 104 °F (40 °C)

Normal operating environment Intended for indoor use only

Maximum altitude 2,000 m (6,562 ft)

www.agilent.com/chem

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2015 Printed in the USA August 18, 2015 5991-6379EN

